SILADENT

Dr. Böhme & Schöps GmbH

Das System

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Instructions for use

Name and general description

Dental blanks made of Zirkon BioStar, Zirkon BioStar Z or Zirkon BioStar Colour are dental materials (semi-finished products) made of yttrium stabilized, pre-sintered zirconium dioxide for milled production of crowns and bridge frameworks on commercial CAD/CAM systems or hand-operated copy-milling machines with outstanding biocompatibility and high resistance against tension and pressure.

The pre-sintered blanks are eminently suitable for machining and have excellent edge stability. On account of the special production process, extremely constant firing shrinkage is achieved, even over different batches, meaning that in most cases it is not necessary to adjust the enlargement factor on the CAD/CAM unit. Only diamond-charged grinding tools or hard metal tools are to be used for machining.

Geometric variants

The blanks are manufactured in various product shapes and sizes. Discs 90 to 110 mm in diameter and 10 to 28 mm in height (as full disc or with flange according to customer requirements). Cuboids and cylinders in various dimensions Other dimensions according to customer requirements are possible by agreement

Material types

The material types Zirkon BioStar white opaque and colour, Zirkon BioStar Z white translucent, Zirkon BioStar K, Zirkon BioStar S differ regarding the strength value achievable after the final sintering and contain a different amount of aluminum oxide. They differ in the manner of whiteness and translucency and both show a different resistance against hydrothermal aging. Zirkon BioStar Colour is already persistent coloured in the pre-sintered stage. Consequently, the dental surgeon or dental technician can select the ideal base material for the respective case. On account of the high strength values, both materials are also suitable for multiple-element posterior bridge-work.

If required, additional physical and chemical key values are referred in our material data sheets. All material types were tested by an accredited testing laboratory regarding biological compatibility and classified as eminently suitable for use in the dental sector.


Indication

Crowns in the anterior and posterior dental area and also bridge frameworks with a maximum of two pontics. When undertaking the preparatory work, it must be ensured that there is an adequate connector surface between the abutment crown and the pontic and adequate minimum wall strengths.

Contraindications

In the case of bruxism, inadequate hard tooth structure or insufficient space, the use of full ceramic crowns and bridges is not recommended or only advised in some cases. The same applies in the case of inadequate oral hygiene or inadequate coordination of masticatory movements.

A full ceramic crown or bridge may also not be used in the event of unsatisfactory preparation results or in the case of production faults.

Ceramic micrograph of zirconium oxide after sintering

Preparation information

The minimum wall strength values in the following table are merely reference values. The decision on the manner and undertaking of the preparation is the responsibility of the dental surgeon and dental technician.

Minimum wall thicknesses:

Double crown primary components	Incisal/occlusal 0.7 mm	Circular 0.5 mm
Single crown framework	Incisal/occlusal 0.7 mm	Circular 0.5 mm
Posterior teeth	Incisal/occlusal 0.7 mm	Circular 0.5 mm
Anterior teeth	Incisal/occlusal 0.7 mm	Circular 0.5 mm
Abutment crowns for bridge frameworks	with 1 pontic	
	Incisal/occlusal 0.7 mm	Circular 0.7 mm
	with 2 pontics	
	Incisal/occlusal 1.0 mm	Circular 0.7 mm

Minimum connector areas:

Anterior bridge framework	with one pontic	7 mm ²
Anterior bridge framework	with two pontics	9 mm ²
Posterior bridge framework	with one pontic	9 mm ²
Posterior bridge framework	with two pontics	12 mm ²
Cantilever bridge framework		12 mm ²

Milling

Milling can be undertaken with most commercially available CAD/CAM milling machines or hand-operated copy-milling machines. The machine manufacturer's manual should be consulted regarding the selection of the milling cutters and the necessary revolutions/minute. In order to avoid staining or a loss of strength, the blank should if possible not come into contact with liquids (water, oils) or fats (handcream) during milling and related handling of the blank.

Wet machining of the dental blanks in their pre-sintered state is in principle possible, but does require that the restoration be dried completely before sintering. Drying can be undertaken in a drying chamber at approx. 80°C or under an infrared lamp. Alternatively, a low-temperature firing program can be selected which will not damage the blank.

Notches, sharp edges or slight damage, also on the underside of the framework, may reduce the strength of the completed sintered work. Consequently, it is already important to ensure that surfaces are well rounded and smooth when milling.

The framework must be optically checked after milling to ensure there are:

- **No discolourations - No flaking - No cracks -**

If any faults are discovered, the framework must not be used to fabricate a restoration.

Colouring

If required, the frameworks may be matched to the desired tooth colour with the standard colouring solutions which are licensed for zirconium oxide ceramics. In order to achieve even results, it must be ensured that the framework is completely clean, fat-free and dry before colouring. The instructions of the manufacturers of the colouring solutions should be followed for the further work procedures.

Sintering

The sintering of the preparations can be undertaken in all standard furnaces which are suitable for zirconium oxide ceramics. The manufacturer's operating instructions must be followed in all cases.

When sintered, the frameworks shrink by around 20 %. The required movement is only possible if the individual crowns or frameworks are fired on corresponding sinter bearings. The framework must be positioned in such a way that it cannot come into contact with other frameworks or the sinter carrier in order to avoid any deformation. The sinter carriers must be examined at least every two weeks in accordance with the manufacturer's instructions to ensure that they are level. Anterior and posterior dental frameworks with three or more elements should generally be positioned with the occlusal side upwards. If sintering wires or pegs are used, it must be ensured that there is no contact with the cap wall area. Alternatively, sintering may be carried out in a sintering crucible in which the restoration is embedded in suitable granular material. One suitable type of granules are ZrO₂ sintering spheres. When using sintering crucibles and granules, the framework is placed in the granules with the occlusal side downwards.

The shrinkage of the preparations made of Zirkon BioStar and Zirkon BioStar Z remains at an extremely constant level, both within the individual batches and over several batches, on account of the special production process. We shall advise you of the exact values with each delivery.

Please refer to the information/parameters in the instructions for use issued by the manufacturer of the equipment you are using.

Depending on the material, the recommended sintering temperatures are:

Sintering program for Zirkon BioStar, white opaque and colour						
	single crown		bridge of 3 - 4 pieces		bridge of > 4 pieces & circular bridges	
	°C	min.	°C	min.	°C	min.
Preheating rate up to 900 °C °C / min.	15	59	12	73	8	110
Holding time at 900 °C	900	5	900	15	900	30
Preheating rate to final temperature of 1.450 °C	15	36	12	46	8	69
Holding time at final temperature of 1.450 °C	1450	30	1450	120	1450	120
Cooling rate down to 900 °C °C / min.	10	55	5	110	4	138
Holding time at 900 °C	900	5	900	10	900	20
Cooling rate down to 200 °C	10	70	5	135	4	175
Total time in h	4 h 20 min.		8 h 29 min.		11 h 02 min.	

Sintering program for Zirkon BioStar Z white translucent						
	single crown		bridge of 3 - 4 pieces		bridge of > 4 pieces & circular bridges	
	°C	min.	°C	min.	°C	min.
Preheating rate up to 900 °C °C / Min.	15	59	12	73	8	110
Holding time at 900 °C	900	5	900	15	900	30
Preheating rate to final temperature of 1.450 °C	15	40	12	50	8	75
Holding time at final temperature of 1.450 °C	1500	30	1500	120	1500	120
Cooling rate down to 900 °C °C / min.	10	60	5	120	4	150
Holding time at 900 °C	900	5	900	10	900	20
Cooling rate down to 200 °C	10	70	5	135	4	175
Total time in h	4 h 29 min.		8 h 43 min.		11 h 20 min.	

Depending on the size and complexity of the individual restoration, firing programs with quicker heating rates may also be selected. Some dental sintering furnaces have fixed firing curves which cannot be altered. The user must check whether the above heating rates are achievable with the furnace being used and are tolerated by the specific restorations without any distortion or crack formation, and define individual settings. As the temperature measurement in every furnace can give somewhat different results, it may possibly be necessary in an individual case to adjust the specific furnace temperature by carrying out a test firing cycle. Massive drops in temperature when opening the sintering furnace should be avoided, as opening the heated furnace can have a negative effect on both the furnace and the sintered products, damaging the restoration. The temperatures vary according to furnace type. With our zirconia products, avoid opening the furnace at temperatures above 200°C. Please refer to the instructions provided by the furnace manufacturer. Dental frameworks and crowns made of zirconium oxide should be fired in a furnace which is only used for these products. If additionally restorations made of other ceramic materials are fired or if liquid ceramics are baked in the same furnace, this may impede the sintering behavior or result in localized discolouration. If the firing of different materials in the same furnace cannot be avoided, then an empty firing cycle at maximum temperature is recommended before firing zirconium oxide ceramics.

Facing of the frameworks

Crowns and bridge frameworks made of Zirkon BioStar materials may be faced with standard facing ceramics for frameworks made of zirconium oxide. The facings passed by the respective manufacturers for zirconium oxide are adapted regarding the thermal expansion coefficient and the material bonding.

Please consult the descriptions of the manufacturers of the facing material for information regarding whether thermal or chemical pretreatment of the sintered construction is recommended before facing.

Finishing after final sintering

If any finishing of the sintered structures is required, they should always be cooled using water and very little pressure applied during finishing to avoid the zirconium oxide undergoing a phase change from tetragonal to monoclinic. It is the tetragonal crystal structure that gives zirconium oxide its high strength. Localized working of the zirconium oxide can create so much pressure or temperature peaks that a phase change takes place. Zirconium oxide in the monoclinic phase has a different density to that in the tetragonal phase and thus also different refraction properties and shade effects.

Information on storage and shelf life

There is no empirical evidence on the life span of presintered dental blanks. However, no change in their material properties is anticipated if they are stored properly, even for long periods. The blanks should not be exposed to direct sunlight. A maximum shelf life therefore cannot and need not be defined. However, as a general precaution, the dental blanks should be processed and the final sintering completed within three years of manufacture.

The minimum wall thicknesses for crown and bridge frameworks and sintering program heating rates given in the instructions for use have been tested many times in practice and have proven satisfactory. However, it may be necessary due to unusual framework size or unequal material concentrations to adjust the processing conditions to the specific restoration. Therefore simply adhering to the listed parameters cannot guarantee accurate results.

Disposal of material residues

Remains of the blanks and, of course, also sintered parts made of Zirkon BioStar materials are categorized as chemically completely inert and are not subject to any special disposal stipulations.

Secondary effects

Such as allergies to zirconium dioxide may very rarely occur.

Reciprocal actions

unknown

Reactions

In case of known incompatibilities and allergies to contents of the zirconium dioxide.

CE-Marking regarding to MDD

The conformity to directive 93/42/EEC is controlled in an assessment procedure by DEKRA Certification GmbH.

Material data sheets for the different Zirkon BioStar materials:

Zirkon BioStar Z

Chemical formula:	ZrO ₂ / Y ₂ O ₃
Chemical name:	Yttrium-stabilized tetragonal zirconium oxide
Colour:	White, translucent
Product description:	Presintered zirconium oxide components for dental use

Chemical composition

ZrO ₂ / HfO ₂ / Y ₂ O ₃	> 99.0 wt%
Y ₂ O ₃	5.15 +/- 0.20 wt%
HfO ₂	< 5.0 wt%
Al ₂ O ₃	< 0.1 wt%
Fe ₂ O ₃	< 0.1 wt %
Na ₂ O	< 0.04 wt %

Physical material data after final sintering

Density	> 6,05 g/cm ³
Bending strength	1.000 +/- 200 MPa
Weibull modulus	> 5
E-modulus	210 GPa
Thermal expansion	11 10 ⁻⁶ 1/K ¹
Open porosity	None

Radioactivity

Activity ²³⁸ U, ²²⁶ Ra, ²³² Th	< 0,2 Bq/g	Meets ISO 13356.2
---	------------	-------------------

Biocompatibility and chemical solubility

Cytotoxicity	Non-toxic DIN EN ISO 10993
Chemical solubility	Below detection limit DIN EN ISO 10903 / DIN EN ISO 6872 (< 1 µg/cm ²)

The above data are typical values and relate only to the test pieces from which they were obtained. They are therefore not unconditionally applicable to other forms of the same material. The values are therefore to be taken as indicative only.

Zirkon BioStar (white opaque) Zirkon BioStar S Zirkon BioStar K

Chemical formula:	ZrO ₂ / Y ₂ O ₃ / Al ₂ O ₃
Chemical name:	Yttrium stabilized tetragonal zirconium oxide
Colour:	White opaque
Product description:	Presintered zirconium oxide components for dental use

Chemical composition

ZrO ₂ / HfO ₂ / Y ₂ O ₃	> 99.0 wt%
Y ₂ O ₃	5.15 +/- 0.20 wt%
HfO ₂	< 5.0 wt%
Al ₂ O ₃	0.25 +/- 0.10 wt%
Fe ₂ O ₃	< 0.1 wt %
Na ₂ O	< 0.04 wt %

Physical material data after final sintering

Density	> 6,05 g/cm ³
Bending strength	1.200 +/- 200 MPa
Weibull modulus	> 8
E-modulus	210 GPa
Thermal expansion	11 10 ⁻⁶ 1/K ¹
Open porosity	None

Radioactivity

Activity ²³⁸ U, ²²⁶ Ra, ²³² Th	< 0,2 Bq/g	Meets ISO 13356.2
---	------------	-------------------

Biocompatibility and chemical solubility

Cytotoxicity	Non-toxic	DIN EN ISO 10993
Chemical solubility	Below detection limit	DIN EN ISO 10903 / DIN EN ISO 6872 (< 1 µg/cm ²)

The above data are typical values and relate only to the test pieces from which they were obtained. They are therefore not unconditionally applicable to other forms of the same material. The values are therefore to be taken as indicative only.

Zirkon BioStar colour 124

Chemical formula:	ZrO ₂ / Y ₂ O ₃
Chemical name:	Yttrium-stabilized tetragonal zirconium oxide
Colour:	Tooth coloured in various shades
Product description:	Presintered zirconium oxide blanks for dental use

Chemical composition

ZrO ₂ / HfO ₂ / Y ₂ O ₃	> 96,0 wt%
Y ₂ O ₃	4,9 - 5,3 wt%
HfO ₂	0 - 5,0 wt%
Al ₂ O ₃	0 - 0,3 wt%
Fe ₂ O ₃	0 - 2,0 wt %
NA ₂ O	0 - 0,04 wt%

Physical material data after final sintering

Density	> 6,05 g/cm ³
Bending strength	1.000 +/- 200 MPa
Weibull modulus	> 5
E-modulus	210 GPa
Thermal expansion	10,5 – 10,8 10 ⁻⁶ 1/K ¹
Open porosity	None

Radioactivity

Activity ²³⁸ U, ²²⁶ Ra, ²³² Th	< 0,2 Bq/g	Meets ISO 13356.2
---	------------	-------------------

Biocompatibility and chemical solubility

Cytotoxicity	Non-toxic
Chemical solubility	Below detection limit DIN EN ISO 10903 / DIN EN ISO 6872 (< 1 µg/cm ²)

The above data are typical values and relate only to the test pieces from which they were obtained. They are therefore not unconditionally applicable to other forms of the same material. The values are therefore to be taken as indicative only.

Further questions: Contact SILADENT technical services (Tel.: +49 (0) 53 21-37 79 25/26) or our sales representatives.

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La valutazione della conformità viene verificata da DEKRA Certification GmbH.